




**SURTEX**  
INSTRUMENTS

*Crafted to Excellence*



AN

## ANESTHESIA

Surtext Laryngoscopes are designed and manufactured to make intubation simple, convenient and as safe as possible. The outcome is a broad product range featuring innovative design, performance, quality and overall value.

Our range features all popular Macintosh & Miller patterns. Our **Corona**™ premium laryngoscopes and **MAXBright**™ LED sets the industry standard for effortless atraumatic intubation and features distal fibre optic (F.O) for bright illumination and maximum vision of the epiglottis and vocal cords. Our **pollo**™ standard features direct bulb illumination utilizing Surtex xenon illuminated bulbs.

Narkose

Anesthesia


Anesthésie

Anesthesia



Catheter Introducing Forceps, Laryngoscopes and Accessories  
 Kathetereinführzangen, Laryngoskope und Zubehör  
 Pinzas para catéteres, Laringoscopios y accesorios  
 Pinces pour introduire les sondes, Laryngoscopes et accessoires  
 Pinze per introdurre le sonde, Laringoscopi ed accessori

**Magill Catheter Introducing Forcep**  
 For Babies



Catalog #	Size
AN-010-15	14.5 cm (5 ¾"), For Babies
AN-010-20	20 cm (8"), For Children
AN-010-25	25 cm (9 ¾"), For Adults

**Apollo™ Standard Miller-Baby Laryngoscope Set**  
 With Battery Handle Ref: AN-290-01 and 2 Blades Ref: AN-110-00 to AN-110-01



Catalog #	Size
AN-100-00	Standard

**Apollo™ Standard Miller Laryngoscope Set**  
 With Battery Handle Ref: AN-290-01 and 3 Blades Ref: AN-110-01 to AN-110-03



Catalog #	Size
AN-101-00	Standard

**Apollo™ Standard Miller Laryngoscope Set**  
 With Battery Handle Ref: AN-290-01 and 4 Blades Ref: AN-110-01 to AN-110-04



Catalog #	Size
AN-102-00	Standard

**Apollo™ Standard Miller Laryngoscope Set**  
 With Battery Handle Ref: AN-290-01 and 5 Blades Ref: AN-110-00 to AN-110-04



Catalog #	Size
AN-103-00	Standard



**Apollo™ Standard Miller Laryngoscope Blade**



Catalog #	Blade Length, Fig, Usage
AN-110-00	55 mm, Fig. 0, For Babies
AN-110-01	80 mm, Fig. 1, For Children
AN-110-02	130 mm, Fig. 2, Adolescents
AN-110-03	170 mm, Fig. 3, Women
AN-110-04	180 mm, Fig. 4, Men



**Apollo™ Standard McIntosh-Baby Laryngoscope Set**

With Battery Handle Ref: AN-290-01 and 2 Blades Ref: AN-210-00 to AN-210-01



Catalog #	Size
AN-200-00	Standard



**Apollo™ Standard McIntosh Laryngoscope Set**

With Battery Handle Ref:- AN-290-01 and 3 Blades Ref:- AN-210-01 to AN-210-03



Catalog #	Size
AN-201-00	Standard



**Apollo™ Standard McIntosh Laryngoscope Set**

With Battery Handle Ref:- AN-290-01 and 4 Blades Ref:- AN-210-01 to AN-210-04



Catalog #	Size
AN-202-00	Standard



**Apollo™ Standard McIntosh Laryngoscope Set**

With Battery Handle Ref:- AN-290-01 and 5 Blades Ref:- AN-210-00 to AN-210-04



Catalog #	Size
AN-203-00	Standard





Laryngoscopes and Accessories  
Laryngoskope und Zubehör  
Laringoscopios y accesorios  
Laryngoscopes et accessoires  
Laringoscopi ed accessori

**Apollo™ Standard McIntosh Laryngoscope Blade**  
Fig. 0, For Babies



Catalog #	Blade Length, Fig, Usage
AN-210-00	55 mm, Fig. 0, For Babies
AN-210-01	75 mm, Fig. 1, For Children
AN-210-02	90 mm, Fig. 2, For Adolescents
AN-210-03	110 mm, Fig. 3, For Women
AN-210-04	135 mm, Fig. 4, For Men



**Apollo™ Standard Battery Handle**  
Chrome Plated



Catalog #	Size
AN-290-01	Standard



**Apollo™ Standard Spare Bulb Large**  
For Blades Fig. 2 to 4



Catalog #	Size
AN-294-01	Standard



**Apollo™ Standard Spare Bulb Small**  
For Blades Fig. 0 to 1



Catalog #	Size
AN-294-02	Standard





Laryngoscopes and Accessories  
Laryngoskope und Zubehör  
Laringoscopios y accesorios  
Laryngoscopes et accessoires  
Laringoscopi ed accessori

**Corona™ Premium Fibre Optic McIntosh-Baby Laryngoscope Set**  
With Battery Handle Ref:- AN-590-01 and 2 Blades Ref:- AN-410-00 to AN-410-01



Catalog #	Size
AN-400-00	Standard



**Corona™ Premium Fibre Optic McIntosh Laryngoscope Set**  
With Battery Handle Ref:- AN-590-01 and 3 Blades Ref:- AN-410-01 to AN-410-03



Catalog #	Size
AN-401-00	Standard



**Corona™ Premium Fibre Optic McIntosh Laryngoscope Set**  
With Battery Handle Ref:- AN-590-01 and 4 Blades Ref:- AN-410-01 to AN-410-04



Catalog #	Size
AN-402-00	Standard



**Corona™ Premium Fibre Optic McIntosh Laryngoscope Set**  
With Battery Handle Ref:- AN-590-01 and 5 Blades Ref:- AN-410-00 to AN-410-04



Catalog #	Size
AN-403-00	Standard



**Corona™ Premium Fibre Optic McIntosh Laryngoscope Blade**



Catalog #	Blade Length, Fig, Usage
AN-410-00	55 mm, Fig. 0, For Babies
AN-410-01	70 mm, Fig. 1, For Children
AN-410-02	80 mm, Fig. 2, For Adolescents
AN-410-03	105 mm, Fig. 3, For Women
AN-410-04	130 mm, Fig. 4, For Men





Laryngoscopes and Accessories  
Laryngoskope und Zubehör  
Laringoscopios y accesorios  
Laryngoscopes et accessoires  
Laringoscopi ed accessori

**Corona™ Premium Fibre Optic Miller-Baby Laryngoscope Set**

With Battery Handle Ref:- AN-590-01 and 2 Blades Ref:- AN-510-00 to AN-510-01



Catalog #	Size
AN-500-00	Standard



**Corona™ Premium Fibre Optic Miller Laryngoscope Set**

With Battery Handle Ref:- AN-590-01 and 3 Blades Ref:- AN-510-01 to AN-510-03



Catalog #	Size
AN-501-00	Standard



**Corona™ Premium Fibre Optic Miller Laryngoscope Set**

With Battery Handle Ref:- AN-590-01 and 4 Blades Ref:- AN-510-01 to AN-510-04



Catalog #	Size
AN-502-00	Standard



**Corona™ Premium Fibre Optic Miller Laryngoscope Set**

With Battery Handle Ref:- AN-590-01 and 5 Blades Ref:- AN-510-00 to AN-510-04



Catalog #	Size
AN-503-00	Standard



**Corona™ Premium Fibre Optic Miller Laryngoscope Blade**

Fig. 0, For Babies



Catalog #	Blade Length, Fig, Usage
AN-510-00	55 mm, Fig. 0, For Babies
AN-510-01	80 mm, Fig. 1, For Children
AN-510-02	130 mm, Fig. 2, For Adolescents
AN-510-03	170 mm, Fig. 3, For Women
AN-510-04	180 mm, Fig. 4, For Men






Laryngoscopes and Accessories  
Laryngoskope und Zubehör  
Laringoscopios y accesorios  
Laryngoscopes et accessoires  
Laringoscopi ed accessori

**Corona™ Premium Fibre Optic Battery Handle**  
With 2.5 Volts Xenon Bulb



Catalog #	Size
AN-590-01	Standard

**Corona™ Premium Fibre Optic Battery Handle**  
Slim Profile, With 2.5 Volts Xenon Bulb



Catalog #	Size
AN-592-02	Standard

**Corona™ Premium Fibre Optic Spare Bulb**  
Xenon



Catalog #	Power Rating
AN-594-01	2.5 V



Laryngoscopes and Accessories  
Laryngoskope und Zubehör  
Laringoscopios y accesorios  
Laryngoscopes et accessoires  
Laringoscopi ed accessori

**MaxBright™ LED Fibre Optic McIntosh-Baby Laryngoscope Set**

With Battery Handle Ref:- AN-890-01 and 2 Blades Ref:- AN-710-00 to AN-710-01



Catalog #	Size
AN-700-00	Standard



**MaxBright™ LED Fibre Optic McIntosh Laryngoscope Set**

With Battery Handle Ref:- AN-890-01 and 3 Blades Ref:- AN-710-01 to AN-710-03



Catalog #	Size
AN-701-00	Standard



**MaxBright™ LED Fibre Optic McIntosh Laryngoscope Set**

With Battery Handle Ref:- AN-890-01 and 4 Blades Ref:- AN-710-01 to AN-710-04



Catalog #	Size
AN-702-00	Standard



**MaxBright™ LED Fibre Optic McIntosh Laryngoscope Set**

With Battery Handle Ref:- AN-890-01 and 5 Blades Ref:- AN-710-00 to AN-710-04



Catalog #	Size
AN-703-00	Standard



**MaxBright™ LED Fibre Optic McIntosh Laryngoscope Blade**

Fig. 0, For Babies



Catalog #	Blade Length, Fig, Usage
AN-710-00	55 mm, Fig. 0, For Babies
AN-710-01	70 mm, Fig. 1, For Children
AN-710-02	80 mm, Fig. 2, For Adolescents
AN-710-03	105 mm, Fig. 3, For Women
AN-710-04	130 mm, Fig. 4, For Men







Laryngoscopes and Accessories  
Laryngoskope und Zubehör  
Laringoscopios y accesorios  
Laryngoscopes et accessoires  
Laringoscopi ed accessori

**MaxBright™ Fibre Optic Miller-Baby Laryngoscope Set**

With Battery Handle Ref:- AN-890-01 and 2 Blades Ref:- AN-810-00 to AN-810-01



Catalog #	Size
AN-800-00	Standard



**MaxBright™ Fibre Optic Miller Laryngoscope Set**

With Battery Handle Ref:- AN-890-01 and 3 Blades Ref:- AN-810-01 to AN-810-03



Catalog #	Size
AN-801-00	Standard



**MaxBright™ Fibre Optic Miller Laryngoscope Set**

With Battery Handle Ref:- AN-890-01 and 4 Blades Ref:- AN-810-01 to AN-810-04



Catalog #	Size
AN-802-00	Standard



**MaxBright™ Fibre Optic Miller Laryngoscope Set**

With Battery Handle Ref:- AN-890-01 and 5 Blades Ref:- AN-810-00 to AN-810-04



Catalog #	Size
AN-803-00	Standard



**MaxBright™ Fibre Optic Miller Laryngoscope Blade**

Fig. 0, For Babies



Catalog #	Blade Length, Fig, Usage
AN-810-00	55 mm, Fig. 0, For Babies
AN-810-01	80 mm, Fig. 1, For Children
AN-810-02	130 mm, Fig. 2, For Adolescents
AN-810-03	170 mm, Fig. 3, For Women
AN-810-04	180 mm, Fig. 4, For Men





Laryngoscopes and Accessories  
Laryngoskope und Zubehör  
Laringoscopios y accesorios  
Laryngoscopes et accessoires  
Laringoscopi ed accessori

**MaxBright™ Fibre Optic Battery Handle**  
With 2.5 Volts LED Bulb




Catalog #	Size
AN-890-01	Standard



**MaxBright™ Fibre Optic Battery Handle**  
Slim Profile, With 2.5 Volts LED Bulb



Catalog #	Size
AN-892-02	Standard



**MaxBright™ LED Fibre Optic Spare Bulb**



Catalog #	Power Rating
AN-894-01	2.5 V

