




**SURTEX**  
INSTRUMENTS

*Crafted to Excellence*



HW

## HOLLOWARE



Surtext holloware products are multipurpose and utilized in medical procedures and hospital environments worldwide. Our holloware is made from the highest-grade non-magnetic 18/8 surgical stainless steel complying with international ISO standards ensuring highly corrosion resistant products that last a lifetime.

All holloware is polished and surface treated to a smooth satin finish to further enhance corrosion resistance and prevent glare and distractions during operative procedures. Our range includes and is not limited to the following bedpans, kidney trays, urinals, surgical instrument trays, gallipot & septum mugs, dressing drums, forceps jars etc.

Hohlwaren

Holloware

Platerie

Holloware



Round Bowls, Flat Dishes  
Runde Schüsseln, Flachschaalen  
Cubetas redondas, Cubetas planas  
Cuvettes rondes, Cuvettes plates  
Scodelle rotonde, Scodelle piatte

**Round Bowl**

Catalog #	Size, Capacity
HW-010-16	Ø 150 x 74 mm, 750 ccm
HW-010-17	Ø 160 x 65 mm, 1000 ccm
HW-010-22	Ø 220 x 70 mm, 2000 ccm
HW-010-24	Ø 240 x 95 mm, 3000 ccm
HW-010-27	Ø 260 x 130 mm, 4500 ccm
HW-010-31	Ø 300 x 105 mm, 6000 ccm
HW-010-34	Ø 320 x 140 mm, 8500 ccm
HW-010-38	Ø 380 x 125 mm, 11000 ccm
HW-010-41	Ø 400 x 150 mm, 14000 ccm



**Round Bowl**

Catalog #	Size, Capacity
HW-020-33	Ø 330 x 90 mm, 6000 ccm
HW-022-33	Ø 330 x 115 mm, 8000 ccm



**Round Bowl**

Catalog #	Size, Capacity
HW-026-11	Ø 110 x 25 mm, 150 ccm
HW-028-11	Ø 110 x 40 mm, 250 ccm



**Lid Only**

Catalog #	Size
HW-030-11	N/A





Stainless Steel Bowls, Measurers  
Metallschalen, Messbecher  
Cubetas de acero inoxidable, Vaso-medidas  
Cuvettes de acier inoxydable, Mesures graduées  
Scodelle di acciaio inossidabile, Bicchieri graduati

**Mixing Bowl/ Solution Bowl**

Catalog #	Size, Capacity
HW-050-04	Ø 40 x 19 mm, 20 ccm
HW-050-06	Ø 61 x 30 mm, 70 ccm
HW-050-08	Ø 80 x 40 mm, 140 ccm
HW-050-11	Ø 116 x 50 mm, 250 ccm
HW-050-12	Ø 128 x 55 mm, 350 ccm
HW-050-14	Ø 147 x 65 mm, 500 ccm
HW-050-16	Ø 167 x 75 mm, 900 ccm
HW-050-18	Ø 187 x 85 mm, 1300 ccm



**Round Dish**

Catalog #	Size, Capacity
HW-060-11	Ø 116 x 35 mm, 200 ccm
HW-060-12	Ø 128 x 38 mm, 300 ccm
HW-060-14	Ø 147 x 40 mm, 400 ccm
HW-060-16	Ø 167 x 45 mm, 500 ccm
HW-060-18	Ø 187 x 45 mm, 800 ccm



**Round Dish**

Catalog #	Size, Capacity
HW-064-11	Ø 110 x 28 mm, 150 ccm



**Graduated Measure**

Catalog #	Capacity
HW-070-01	100 ccm





Measurers, Universal Trays, Instrument Trays  
 Messbecher, Mehrzweckschalen, Instrumententabletts  
 Vaso-medidas, Bandejas universales, Bandejas para instrumentos  
 Mesures graduées, Plateaux universaux, Plateaux pour instruments  
 Bicchieri graduati, Bacinelle universali, Bacinelle per strumenti

**Graduated Measure**

Catalog #	Capacity
HW-072-01	250 ccm
HW-072-02	500 ccm
HW-072-03	1000 ccm
HW-072-04	2000 ccm



**Petri Dish**

Catalog #	Size, Profile
HW-090-00	Ø 100 x 20 mm, With Lid
HW-090-01	Ø 96 x 20 mm, Without Lid
HW-090-02	N/A, Lid Only



**Universal Tray**

Catalog #	Size
HW-100-23	230 x 155 x 30 mm
HW-100-30	300 x 175 x 30 mm
HW-100-31	300 x 175 x 40 mm
HW-100-35	350 x 240 x 36 mm
HW-100-40	470 x 320 x 50 mm



**Instrument Tray**

Catalog #	Size
HW-110-21	210 x 160 x 10 mm
HW-110-24	240 x 180 x 10 mm
HW-110-30	310 x 150 x 10 mm
HW-110-31	310 x 210 x 10 mm
HW-110-35	350 x 240 x 10 mm
HW-110-40	400 x 270 x 20 mm
HW-110-42	420 x 280 x 10 mm





Universal Trays, Instrument Trays, Metal Jars  
Mehrzweckschalen, Instrumententablets, Standgefäße  
Bandejas universales, Bandejas para instrumentos, Tarros metálicos  
Plateaux universaux, Plateaux pour instruments, Récipients en métal  
Bacinelle universali, Bacinelle per strumenti, Recipienti metallici

**Universal Tray**



Catalog #	Size
HW-120-23	235 x 190 x 40 mm

**Instrument Tray  
Without Cover**



Catalog #	Size
HW-122-01	150 x 100 x 25 mm
HW-122-02	200 x 150 x 40 mm
HW-122-03	250 x 200 x 40 mm
HW-122-04	300 x 250 x 50 mm
HW-122-05	350 x 250 x 50 mm
HW-122-06	350 x 300 x 50 mm
HW-122-07	400 x 300 x 50 mm
HW-122-08	450 x 350 x 50 mm

**Instrument Tray  
With Cover**



Catalog #	Size
HW-124-01	150 x 100 x 25 mm
HW-124-02	200 x 150 x 40 mm
HW-124-03	250 x 200 x 40 mm
HW-124-04	300 x 250 x 50 mm
HW-124-05	350 x 250 x 50 mm
HW-124-06	350 x 300 x 50 mm
HW-124-07	400 x 300 x 50 mm
HW-124-08	450 x 350 x 50 mm

**Round Bowl**



Catalog #	Size
HW-130-22	Ø 220 x 70 mm



Instrument Boxes  
Instrumentenkästen  
Cajas para instrumentos  
Boîtes pour instruments  
Scatole per strumenti

**Metal Jar**

Catalog #	Size
HW-150-30	Ø 30 x 90 mm
HW-150-50	Ø 50 x 100 mm
HW-150-51	Ø 50 x 130 mm
HW-150-52	Ø 50 x 175 mm
HW-150-75	Ø 75 x 100 mm
HW-150-76	Ø 75 x 130 mm
HW-150-77	Ø 75 x 175 mm
HW-150-98	Ø 100 x 175 mm
HW-150-99	Ø 100 x 200 mm



**Lid Only**

Catalog #	Diamete, Profile
HW-153-50	Ø 50 mm, Closed
HW-153-75	Ø 75 mm, Closed
HW-153-99	Ø 100 mm, Closed
HW-155-50	Ø 50 mm, With Triangular Opening
HW-155-75	Ø 75 mm, With Triangular Opening
HW-155-99	Ø 100 mm, With Triangular Opening



**Instrument Box  
Lid With Knob**

Catalog #	Size
HW-200-16	165 x 85 x 35 mm
HW-200-20	200 x 100 x 35 mm
HW-200-21	200 x 100 x 50 mm
HW-200-23	230 x 130 x 50 mm
HW-200-26	260 x 150 x 50 mm
HW-200-27	270 x 235 x 50 mm
HW-200-30	300 x 200 x 50 mm
HW-200-33	325 x 275 x 50 mm





Instrument Boxes, Round Metal Boxes, Medicine Cups  
 Instrumentenkästen, Runde Metall Dosen, Medizinbecher  
 Cajas para instrumentos, Cajas metálicas redondas, Vasos medicinales  
 Boîtes pour instruments, Boîtes pour instruments, Gobelets  
 Scatole per strumenti, Scatole per strumenti, Bicchieri per medicina

**Instrument Box**  
Lid Without Knob

Catalog #	Size
HW-220-16	165 x 85 x 35 mm
HW-220-20	200 x 100 x 35 mm
HW-220-21	200 x 100 x 50 mm
HW-220-23	230 x 130 x 50 mm
HW-220-26	260 x 150 x 50 mm
HW-220-27	270 x 235 x 50 mm
HW-220-30	300 x 200 x 50 mm
HW-220-32	324 x 275 x 50 mm
HW-220-33	325 x 275 x 50 mm
HW-220-42	420 x 80 x 40 mm
HW-220-43	420 x 175 x 50 mm



**Instrument Box**

Catalog #	Size
HW-230-20	205 x 105 x 35 mm



**Dressing Jar**

Catalog #	Size
HW-250-06	Ø 75 x 60 mm
HW-250-07	Ø 75 x 75 mm
HW-252-06	Ø 100 x 60 mm
HW-252-10	Ø 100 x 100 mm
HW-254-08	Ø 125 x 80 mm
HW-254-12	Ø 125 x 125 mm
HW-256-15	Ø 150 x 150 mm





Cups  
Becher  
Tarros  
Réipients  
Recipienti

**Medicine Cup**

Catalog #	Capacity
HW-300-25	25 cc
HW-300-50	50 cc



**Medicine Spoon**

Catalog #	Capacity
HW-310-20	20 cc



**Tablet Tray**

Catalog #	Size
HW-320-75	Ø 75 x 12 mm



**Metal Cup  
With Lid**

Catalog #	Size
HW-350-90	Ø 90 x 90 mm



**Denture Cup  
With Lid**

Catalog #	Size
HW-360-75	Ø 75 x 50 mm







Cups, Kidney Dishes, Bed Pans  
Becher, Nierenschale, Stechbecken  
Tarros, Riñoneras, Silletas  
Récipients, Bassins réniformes,  
Recipienti, Bacinelle reniformi, Padelle

**Sputum Cup**

Catalog #	Size
HW-370-75	Ø 75 x 80 mm



**Kidney Dish**

Catalog #	Size
HW-400-17	170 x 100 x 35 mm
HW-400-20	250 x 140 x 38 mm
HW-400-25	250 x 140 x 40 mm
HW-400-27	275 x 150 x 45 mm



**Lid For Kidney Dish**

Catalog #	Size
HW-402-25	252 mm



**Bed Pan**

Catalog #	Size
HW-500-23	Ø 235 mm
HW-500-31	Ø 310 mm



**Bed Pan**

Catalog #	Size
HW-510-00	N/A






Triangular Basins, Placenta Pan, Cotton Dispensers  
 Dreieckschale, Plazentabecken, Wattespende  
 Cubetas triangulares, Cubeta para placenta, Dispensadores de algodón  
 Bassins triangulaires, Bassin pour placenta, Distributeurs de coton  
 Bacinelle triangolari, Bacinella per placenta, Distributori di cotone

**Male Urinal**

Catalog #	Size
HW-512-00	N/A



**Female Urinal**

Catalog #	Size
HW-514-00	N/A




**Graduated Triangular Basin  
With Handle**

Catalog #	Capacity
HW-520-30	600 ml



**Graduated Triangular Basin  
Without Handle**

Catalog #	Capacity
HW-522-30	600 ml



**Placenta Pan**

Catalog #	Size
HW-530-36	N/A





Spirit Lamps  
Spirituslampen  
Lámparas de alcohol  
Lampes à alcool  
Lampade a spirito

**Methot Cotton Dispenser**  
With Lid

Catalog #	Size
HW-610-75	Ø 50 x 75 mm



**Methot Cotton Dispenser**  
Without Lid

Catalog #	Size
HW-612-75	Ø 50 x 75 mm



**Spirit Lamp**

Catalog #	Capacity
HW-700-00	120 ccm



**Spirit Lamp**

Catalog #	Capacity
HW-702-00	120 ccm



**Spirit Lamp**

Catalog #	Capacity
HW-704-00	120 ccm





