




SURTEX
INSTRUMENTS

Crafted to Excellence



SF



COTTON SWAB AND DRESSING FORCEPS

Surtex sponge forceps are the perfect choice for holding swabs or sponges in medical procedures. Our range is manufactured in box-joint scissor action with ratcheted handles and looped, smooth or serrated jaws, curved or straight to meet your specific purpose. Many surgeons prefer to use Surtex sponge forceps because they help prevent trauma during surgery by enabling the physician to hold and move tissue very precisely and gently.

Tupfer und Kornzangen

Pizas para tapones y -curaciones

Pinces à tampons et à pansements

Pinze per tamponi e per medicazione



Dressing Forceps
Tupferzangen
Pinzas para tapones
Pinces à tampons
Pinze per tamponi

Ulrich Dressing Forcep
Straight

Catalog #	Size
SF-100-22	22.5 cm (8 ¾")
SF-100-27	26.5 cm (10 ½")



Ulrich Dressing Forcep
Curved

Catalog #	Size
SF-101-22	22.5 cm (8 ¾")
SF-101-27	26.5 cm (10 ½")




Bryant Dressing Forcep
Curved

Catalog #	Size
SF-111-13	13 cm (5")



Lister Sinus Forcep
Straight

Catalog #	Size
SF-120-13	13 cm (5")
SF-120-16	16 cm (6 ¼")
SF-120-18	17 cm (7")



Martin Dressing Forcep
Straight

Catalog #	Size
SF-130-17	17.5 cm (7")






Dressing Forceps
Tupferzangen
Pinzas para tapones
Pincas à tampons
Pinze per tamponi

SF

COTTON SWAB AND DRESSING FORCEPS

Martin Dressing Forcep
Curved




Catalog #	Size
SF-131-17	17 cm (6 3/4")

Gross-Maier Dressing Forcep
Straight, With Ratchet



Catalog #	Size
SF-150-20	20.5 cm (8")
SF-150-22	22 cm (8 3/4")
SF-150-27	26.5 cm (10 1/2")

Gross-Maier Dressing Forcep
Curved, With Ratchet



Catalog #	Size
SF-151-20	20 cm (8")
SF-151-22	21.5 cm (8 1/2")
SF-151-27	26 cm (10 1/4")

Gross-Maier Dressing Forcep
Straight, Without Ratchet



Catalog #	Size
SF-152-20	20.5 cm (8")
SF-152-22	22 cm (8 3/4")
SF-152-27	26.5 cm (10 1/2")



Dressing Forceps
Tupferzangen
Pinzas para tapones
Pinces à tampons
Pinze per tamponi


Gross-Maier Dressing Forcep
Curved, Without Ratchet

Catalog #	Size
SF-153-20	20 cm (8")
SF-153-22	21.5 cm (8 1/2")
SF-153-27	26 cm (10 1/4")




Pelkmann Foreign Body Forcep
Straight

Catalog #	Size
SF-160-26	26 cm (10 1/4")



Pelkmann Foreign Body Forcep
Curved

Catalog #	Size
SF-161-26	26 cm (10 1/4")




Duplay Sponge Holding Forcep
Straight

Catalog #	Size
SF-164-21	21 cm (8 1/4")



Duplay Sponge Holding Forcep
Curved

Catalog #	Size
SF-165-21	21 cm (8 1/4")





Dressing Forceps
Tupferzangen
Pinzas para tapones
Pincas à tampons
Pinze per tamponi

Tunneling Forcep
Curved




Catalog #	Size
SF-171-40	39.5 cm (15 1/2")

Collin Dressing Forcep
Straight




Catalog #	Size
SF-180-25	25.5 cm (10")

Cheron Dressing Forcep
Straight




Catalog #	Size
SF-191-25	23.5 cm (9 1/4")

Emmet Dressing Forcep
Straight



Catalog #	Size
SF-196-21	21 cm (8 1/4")

Bozemann Sponge Holding Forcep
Straight




Catalog #	Size
SF-200-26	26 cm (10 1/4")



Dressing Forceps
Tupferzangen
Pinzas para tapones
Pinces à tampons
Pinze per tamponi


Bozemann Sponge Holding Forcep
Curved S Shaped

Catalog #	Size
SF-201-26	26 cm (10 1/4")



Bozemann-Douglas Sponge Holding Forcep
Curved S Shaped

Catalog #	Size
SF-203-26	26 cm (10 1/4")




Heywood-Smith Sponge Holding Forcep
Straight

Catalog #	Size
SF-220-20	20 cm (8")



Heywood-Smith Sponge Holding Forcep
Straight

Catalog #	Size
SF-222-21	21 cm (8 1/4")



Bergmann Sponge Holding Forcep
Straight

Catalog #	Size
SF-250-24	24 cm (9 1/2")





Dressing Forceps
Tupferzangen
Pinzas para tapones
Pinces à tampons
Pinze per tamponi

Foerster Sponge Holding Forcep
Straight, Smooth Jaws

Catalog #	Size
SF-300-25	25 cm (9 3/4")



Foerster Sponge Holding Forcep
Curved, Smooth Jaws

Catalog #	Size
SF-301-25	24.5 cm (9 3/4")




Foerster Sponge Holding Forcep
Straight

Catalog #	Size
SF-350-18	18 cm (7")
SF-350-20	20.5 cm (8")
SF-350-25	25 cm (9 3/4")



Foerster Sponge Holding Forcep
Curved

Catalog #	Size
SF-351-18	17.5 cm (7")
SF-351-20	20 cm (8")
SF-351-25	24.5 cm (9 3/4")






Dressing Forceps
Tupferzangen
Pinzas para tapones
Pinces à tampons
Pinze per tamponi


Foerster Sponge Holding Forcep
Straight

Catalog #	Size
SF-360-25	25 cm (9 3/4")



Foerster Sponge Holding Forcep
Curved

Catalog #	Size
SF-361-25	24.5 cm (9 3/4")



Ramsey Sponge Holding Forcep
Straight

Catalog #	Size
SF-380-18	18.5 cm (7 1/4")
SF-380-25	24.5 cm (9 3/4")

